

Elcam Medical



Injectable Drug Delivery Devices



Improving 
Quality
of Life



Flexi-Q PFS

Fully Disposable Auto-Injector for Pre-Filled Syringes

The Flexi-Q PFS is a fully disposable auto-injector for Pre-Filled Syringes (PFSs). This single-use auto-injector features passive needle shielding, provides life cycle management capabilities for new and existing drugs, encourages better patient compliance and improves safety in self-injection.

Flexi-Q PFS

Fully Disposable Auto-Injector for Pre-Filled Syringes

Pharmaceutical Company Advantages

Differentiation

- Enables using a drug delivery device for drugs in pre-filled syringes
- Flexible delivery volumes of 0.3 - 1.0 ml.
- Reduced skin reactions - injection starts only after full needle penetration is achieved
- Delivery of entire dose - same dead space as with manual techniques
- Color coding flexibility and other customization options

Contact us for an auto-injector for high-viscosity drugs in PFS.

Short Regulatory Path

- Compatible with standard glass syringes - no need to change drug container

Patient Advantages

Ease of Use - Better Patient Compliance


- Fully automatic injection process
- Quiet injection reduces patients' anxiety and perceived pain
- Reduces needle phobia - by hiding the needle through all the injection process

Patient Safety


- Prevents needle-stick injuries - needle is protected through all the injection process
- Reduced local skin reactions - injects only into the correct depth (SC)
- Clear visual and audible confirmation of the injection progress from beginning to end
- Prevents inadvertent activation - allows activation only after pressing the injector against the patient's skin
- Clear distinction between injector ends - prevents injection in wrong direction

Full Needle Protection


Before Injection



During Injection



After Injection



Inject Button

Drug Pre-Filled Syringe

Long Observation Window

Needle Guard



Needle Cap

Product commercial availability is expected during the first half of 2009

Improving
Quality
of Life



Head Office / Manufacturer
Elcam Medical
BarAm 13860, Israel
Tel: 972-4-698-8120/1/2
Cell: 972-52-831-3890
Fax: 972-4-698-0777
sales@elcam.co.il
www.elcam-medical.com

USA
Elcam Medical INC.
2 University Plaza, Suite 620
07601 Hackensack, NJ
Tel: +1 201 4571120
Tel: +1 201 4571125
info@elcam-medical.com
www.elcam-medical.com

Europe
MedNet GmbH
Borkstrasse 10
48163 Muenster Germany
Tel +49 (0) 251 32266 0
Tel +49 (0) 251 32266 22
info@medneteuropa.com
www.medneteuropa.com

China
Elcam China Office
Room 25B, Building B, No. 1118 of Changshou Rd.,
Shanghai 200042, China
上海市长寿路1118号悦达国际大厦B幢25B室, 200042
Tel: +86 2152375561
Fax: +86 2152375562
johnliu02@126.com