

Elcam Medical, the leading OEM provider of stopcocks continues to provide innovation with the introduction of two new 510(k) cleared and CE approved products – TripleSafe and DoubleSafe Swabbable, Antimicrobial Stopcocks.



**DoubleSafe (DSS)** – Silver ions embedded in the Swabbable Stopcock provides TWO layers of protection:

- 1) A closed system designed to prevent contaminants' penetration into the stopcock fluid path and thereby reducing the risk of bloodstream infections.
- 2) Silver ions to reduce and prevent bacteria colonization on the product.

**TripleSafe (TSS)** – Silver ions embedded in the Marvelous stopcock provides THREE layers of protection:

- 1) A closed system
- 2) Silver ions – just like the DSS configuration
- 3) A unique flushing channel to reduce the residual volume and minimize dead space that further contributes to preventing bacterial colonization

Elcam Medical Silver ion based antimicrobial technology was laboratory tested and proven to be effective against the most prevalent sources of Catheter Related Blood Stream Infections (according to the CDC). The tests showed at least 2 log bacteria reduction and in some bacteria such as: Staph aureus, Staph epidrmidis, Enrerobacter and Candida spp. even more than 4 log bacteria reduction was achieved as well as almost 3 log bacteria reduction against MRSA (**Methicillin Resistant *Staphylococcus aureus***) which is known as a most resistant microorganism.

The new CMS reimbursement ruling on hospital-acquired infections went into effect on October 1, 2008. As leaders in innovation and committed to fighting the war against infections, Elcam Medical is proud to introduce its antimicrobial technologies to its Swabbable Stopcocks.

Elcam Medical is the leading worldwide OEM supplier of fluid management, drug delivery, and vital signs monitoring systems and devices.

For further information, please contact us at [sales@elcam.co.il](mailto:sales@elcam.co.il)